

Medipal® Alcohol wipes are intended for the disinfection of non-porous hard surfaces, equipment and non invasive medical devices within healthcare environments.

Medipal® Alcohol wipes are **proven to kill** most microorganisms such as bacteria, fungi, yeasts and viruses. The wipes have been **independently tested** to a range of **EN standards**.

Bactericidal		
Organism	Contact time	Testing standard
<i>Acinetobacter baumannii</i>	30 sec	EN13697
		EN14561
<i>Enterococcus hirae</i>	30 sec	EN13697
		EN14561
	1 min	EN13727
		EN1276
<i>Escherichia coli</i>	30 sec	EN13697
	1 min	EN1276
MRSA	30 sec	EN13697
		EN14561
<i>Pseudomonas aeruginosa</i>	30 sec	EN13697
		EN14561
	1 min	EN13727
		EN1276
<i>Staphylococcus aureus</i>	30 sec	EN13697
		EN14561
	1 min	EN13727
		EN1276
		EN16615
VRE	30 sec	EN13697
		EN14561

Mycobactericidal		
Organism	Contact time	Testing standard
<i>Mycobacterium avium</i>	1 min	EN14348 EN14563
<i>Mycobacterium terrae</i>	1 min	EN14348 EN14563

Yeasticidal		
Organism	*Contact time	Testing standard
<i>Candida albicans</i>	1 min	EN1650
		EN13624
		EN16615

Virucidal		
Organism	Contact time	Testing standard
Murine Norovirus	1 min	EN 14476
Hepatitis C	5 min	BluTest testing protocol
HIV	5 min	BluTest testing protocol

