

## GP Gynaecological Plinth


Featuring improved access for the clinician and enhanced comfort and dignity for the patient. Effortless adjustment is ensured by electric motors for height and tilt. The electric seat also enables Trendelenburg. Our operators stool is a key optional accessory which can be coordinated with the plinth.

### Standard Features:

- 260kg safe working load for electric mechanism
- Excellent lithotomy access and foot-room for the clinician
- Removable leg section allows use as a standard 2 section plinth
- Retractable rubber castors & Adjustable foot for uneven floors
- Stainless steel accessory rails with a pair of indexing clamps
- Stainless steel tray & choice of patient feet/leg rests included
- Gas spring assisted back section
- MRSA resistant anti bacterial & fire retardant upholstery
- World leading Linak 'High Reliability' actuators
- Upholstery available in 14 standard colours
- Manufactured in the UK and CE Marked
- **Warranty:** 7 Years on Couch Frame & 3 Years on moving parts

### Popular Options & Accessories:

- Fold Down Loop Arms (Shown in image - not included as standard)
- See full range of Options & Accessories on page 30/31



### Plinth Mechanism

Electric Plinth with Stirrups

### Product Code

GY02ES

Electric Plinth with Troughs (Shown)

GY02ET

Electric Plinth with Footpads

GY02EF

## Ultrasound Couch


Provides smooth and effortless patient positioning with electric variable height, electric tilt (providing Trendelenburg facility) and gas spring assisted back and leg sections.

### Standard Features:

- 260kg safe working load for electric mechanism
- Electric tilt facility to allow for Trendelenburg positioning
- Retractable rubber castors
- Adjustable foot for uneven floors
- MRSA resistant anti bacterial & fire retardant upholstery
- Gas spring assisted head & foot sections
- Extra low entry height for wheel chair access
- World leading Linak 'High Reliability' actuators
- Upholstery available in 14 standard colours
- Manufactured in the UK and CE Marked
- **Warranty:** 7 Years on Couch Frame & 3 Years on moving parts

### Popular Options & Accessories:

- Breathing Hole
- Electric Backrest
- Adjustable Arm Rests
- Fold Down Padded Arm Rests (Day-Care Style)
- See full range of Options & Accessories on page 30/31



### Plinth Mechanism

Electric Plinth

### Product Code

US03TE