


HOSPITAL CARE

# Terumo Syringe

## 3-part Syringe Product Information



Disposable syringes  
for single use

Every single detail is the result of Terumo's long time experience in this field which is appreciated by users throughout the world.

- Siliconisation for optimal plunger mobility
- The integral gasket ensures a firm and distinct end-stop
- Clear graduated scale allows accurate dosage
- Transparent barrel for visualisation
- Finger grips to support easy handling

Description	Raw Materials
1 Barrel	Polypropylene
2 Plunger	Polypropylene
3 Gasket	Thermoplastic elastomer (for ,SS'-codes) Polyisoprene rubber (for ,MDSS'-codes)
Lubricant	Silicone oil

Does not contain components made of PVC or natural rubber latex.

### Sterilization

With radiation (e-beam) (for ,SS'-codes) or ethylene oxide (for ,MDSS'-codes) to a Sterility Assurance Level of at least 10<sup>-6</sup>

### Storage

- Shelf life of 5 years
- The batch number and expiry date are printed on the paper (for ,MDSS'-codes) or embossed in the film of the blister pack (for ,SS'-codes)
- Do not store at extreme temperature and humidity (for ,MDSS'- and ,SS\*'-codes) or store between 2°C and 30°C (for ,SS+ '-codes)

## Product Range

Product Code	Capacity		Increments	Tip	Replacement for
	Nominal	Graduated			
MDSS01SE	1 ml	1 ml	0.01 ml	Concentric Luer Slip	SS+01T1
MDSS01LE	1 ml	1 ml	0.01 ml	Concentric Luer Lock	-
SS*02LE1	2.5 ml	2.5 ml	0.1 ml	Concentric Luer Lock	-
MDSS03SE	3 ml	3 ml	0.1 ml	Concentric Luer Slip	SS*02SE1, SS+03S
MDSS03LE	3 ml	3 ml	0.1 ml	Concentric Luer Lock	SS*03L, SS+03L1
MDSS05LE <sup>1</sup>	5 ml	5.5 ml	0.2 ml	Concentric Luer Lock	-
SS*05LE1	5 ml	6 ml	0.2 ml	Concentric Luer Lock	-
MDSS05SE	5 ml	5.5 ml	0.2 ml	Concentric Luer Slip	SS*05SE1, SS+05S
MDSS10SE	10 ml	11 ml	0.2 ml	Concentric Luer Slip	-
MDSS10ESE	10 ml	11 ml	0.2 ml	Eccentric Luer Slip	SS+10ES1
MDSS10LE <sup>1</sup>	10 ml	11 ml	0.2 ml	Concentric Luer Lock	-
SS*10LE1	10 ml	12 ml	0.2 ml	Concentric Luer Lock	-
MDSS20ESE	20 ml	22 ml	1 ml	Eccentric Luer Slip	SS+20ES1
MDSS20LE <sup>1</sup>	20 ml	22 ml	1 ml	Concentric Luer Lock	-
SS+20L1	20 ml	25 ml	1 ml	Concentric Luer Lock	-
SS+30ES1	30 ml	35 ml	1 ml	Eccentric Luer Slip	-
SS+30L1MP	30 ml	35 ml	1 ml	Concentric Luer Lock	-
SS+50L1	50 ml	60 ml	1 ml	Concentric Luer Lock	-
SS+50ES1	50 ml	60 ml	1 ml	Eccentric Luer Slip	-
SS+50C1	50 ml	60 ml	1 ml	Catheter Tip	-



STERILE EO  
for ,MDSS'-codes

STERILE R  
for ,SS'-codes

<sup>1</sup>This concentric luer lock reference is not suited for use with syringe pumps.

## Precautions

- Do not use if package is damaged
- Do not reuse
- Dispose of sharps after single use according to facility protocol to avoid risk of infection

Always consult the product label for a complete overview of warnings, cautions and/or precautions prior to actual use.

Choose from a range of volumes from 1 ml to 50 ml

Syringes are available in all common tip types.

## Packaging

Unit Box	Shipping Carton		
	Quantity	Gross Weight	Dimensions (L x W x H)
100	3000	12.70 kg	635 x 295 x 505 mm
100	3000	13.00 kg	655 x 300 x 520 mm
100	2400	10.40 kg	397 x 581 x 444 mm
100	1800	10.30 kg	445 x 447 x 449 mm
100	1800	11.10 kg	445 x 447 x 449 mm
100	1800	13.20 kg	655 x 418 x 426 mm
100	1600	10.60 kg	412 x 573 x 444 mm
100	1800	12.80 kg	656 x 419 x 426 mm
100	1200	13.00 kg	569 x 473 x 459 mm
100	1200	14.20 kg	568 x 496 x 460 mm
100	1200	13.00 kg	575 x 455 x 455 mm
100	800	15.00 kg	359 x 581 x 383 mm
50	400	8.30 kg	468 x 390 x 455 mm
50	400	8.10 kg	465 x 380 x 445 mm
50	600	9.90 kg	468 x 390 x 455 mm
25	200	6.30 kg	560 x 398 x 255 mm
25	200	6.20 kg	558 x 395 x 254 mm
25	100	5.10 kg	464 x 320 x 275 mm
25	100	5.10 kg	472 x 324 x 282 mm
25	100	5.30 kg	472 x 324 x 282 mm

Single unit: Blister pack  
(film/film for ,SS'-codes and film/paper for ,MDSS'-codes)

## Regulatory Status

CE Certification –  
for ,SS'-codes Class IIa medical device  
for ,MDSS'-codes Class Is medical device

**Notified body**  
TÜV Rheinland 0197

**Manufacturer**  
for ,SS+/-codes TERUMO (Philippines)  
CORPORATION, 124 East Main  
Avenue, Laguna Technopark, Binan,  
Philippines  
for ,SS+/-codes TERUMO CORPORATION,  
2-44-1 Hatagaya, Shibuya-Ku,  
Toyko 151-0072, Japan  
for ,MDSS'-code Yangzhou Medline  
Industry Co., Ltd. (China)

**EC representative**  
for ,SS'-codes TERUMO EUROPE N.V. –  
Leuven, Belgium  
for ,MDSS'-code Shanghai International  
Corp. GmbH (Europe) – Hamburg,  
Germany

**Distributor**  
for ,MDSS'-codes TERUMO EUROPE  
N.V., Interleuvenlaan 40, 3001 Leuven,  
Belgium

Your Terumo partner



**TERUMO EUROPE N.V.**

Terumo Corporation Terumo Europe N.V.  
+81 3 6742 8500 +32 16 38 12 11

www.terumo-europe.com

All brand names are trademarks or registered trademarks owned by TERUMO CORPORATION, its affiliates, or unrelated third parties.  
Published by Terumo Europe N.V.

This brochure is only applicable for countries within EMEA, excluding the MENA region.