



Solid Door



Glass Door

Upright Control Plus Solid/Glass Door Refrigerator 273L

GDPA code: PPSR273UK (Solid), PPGR273UK (Glass)

SKU code: 444410651 (Solid), 444410652 (Glass)

- ✓ 2°C to 8°C temperature range
- ✓ Air and core temperature display
- ✓ Represents your vaccine vial
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Audible and visual temperature alarms
- ✓ 5 year warranty as standard



Specification: Upright Control Plus Refrigerator 273L

Gross Capacity (litres)	273L
Net Capacity (litres)	269L
Unit External Width (mm)	595mm
Unit External Depth (mm)	625mm
Unit External Height (mm)	1565mm
Product Weight (kg)	60kg (solid), 68.5kg (glass)
Packed Width (mm)	650mm
Packed Depth (mm)	727mm
Packed Height (mm)	1585mm
Packed Weight (kg)	65kg (solid), 73.5kg (glass)
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	480mm W x 1382mm H x 445mm
Temperature Range	+2°C to +8°C (tolerance + or - 2)
Colour	White/Green
External Finish	Painted Steel (solid), Painted aluminium/glass (glass)
Internal Finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS mix)
Door Finish	Painted Steel / Solid
Door Type	Solid (solid), Painted aluminium/glass (glass)
Door Style	Hinged, Reversible (solid), Hinged, right handed only, non-reversible (glass)
Door Open Depth (mm)	1185mm
Replacement Door Seal	No (solid), Yes (glass)
Open Door Sensor	Yes
Lock	Yes, 2 Keys
Adjustable Feet (number & position)	2 adjustable front feet, 2 rear rollers
Porthole	No
Porthole Size	N/A
Porthole Location	N/A
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Remote Alarm Contacts	Yes
Audible Alarms	Yes
Visual Alarms	Yes
Number of Shelves	5
Number of Shelf Positions	12
Shelf Size (mm) (WxDxH)	4 x 385mm W x 6mm H x 500mm D 1 x 340mm x 6mm H x 497mm D
Shelf Material	Wire
Number of Baskets	1
Basket Material	Wire
Basket Size (mm) (WxDxH)	—
Interior Light	Yes
Light Type	Halogen
Light Activation	Controller
Fan	Yes
Spark Free	No
Auto or Manual Defrost	Auto
Controller Type	Digital, External
Number of Probes	2
Probe Type	Air & Load (Glycol Solution)
Climate Class	SN
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	Yes

Specification: Upright Control Plus Refrigerator 273L

HCFC Free	Yes
Ammonia Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Noise Level (dB)	41 (solid), 42 (glass)
Energy Consumption (kwh)	0.93kW/24h (solid), 1.08kW/24h (glass)
Voltage Supply	220-240
Aux Output	Yes
Insulation Thickness (mm)	60mm
Insulation Type	Cyclopentane
Blowing Agent	Yes
Anti Microbial Handle	No
Optional Drawer System Available	No (solid), Yes (glass)
Wall Brackets	No
Digital Lock	No