

# Pocket BVM, PBVM-G

Collapsible Bag Valve Mask • Green Case

- Full 1.6L adult bag valve mask
- Contoured adult mask (size 4) provides a tight seal
- Delivers 500-600cc
- Includes 6.5ft of O<sub>2</sub> tubing
- Packed size reduction of 75% compared to standard BVMs
- Can be stored unconventionally
- PEEP compatible



## Delivering Resuscitation with Reduced Size and Weight

The Pocket BVM is a disposable, collapsible, medical-grade silicone resuscitator for adult patients. The device expands from a disc-like hard shell case to a full-size BVM.

The Pocket BVM reduces cube space by up to 75% over traditional BVMs. It is easily stored in unconventional locations and has a ruggedized case capable of withstanding approximately 160lbs psi of direct force. The bag inflates quickly, the mask accommodates most adults, and the unit includes 6.5ft of oxygen tubing. The integrated mask utilizes a curved mold design that mimics the feel and function of inflatable masks without the potential for deflation and delivers a proper seal with each use.

# Pocket BVM, PBVM-G

## Product Highlights

Full-size 1.6L bag

Folds compactly into a hard shell case

Positive End-Expiratory Pressure (PEEP)  
compatible



### This compact, full-size bag valve mask is an ideal component of any medic's kit

- Full-size 1.6L adult bag valve mask (standard 15/22 connection)
- Positive End-Expiratory Pressure (PEEP) compatible
- Delivers 500-600cc
- Size and weight reduction of 75% compared to traditional BVMs
- Can be stored unconventionally compared to other BVMs
- Contoured silicone adult size 4 mask
- Transparent mask provides a clear view of the procedure
- Hard shell case withstands up to 160lbs of direct force
- 6.5ft of robust O<sub>2</sub> tubing included

## Specifications

<b>Description</b> .....	Pocket bag valve mask, 6.5ft O <sub>2</sub> tubing, green case
<b>Part Number</b> .....	PBVM-G
<b>NSN</b> .....	6515-01-593-4841
<b>Packed Size</b> .....	2.75in H x 5.25in D (7cm H x 13.3cm D)
<b>Weight</b> .....	17oz (482g)
<b>Resuscitator Volume</b> .....	1600ml
<b>Body Mass Range</b> .....	>40kg (88lbs)
<b>Patient Connector</b> .....	Inside 15mm female, outside 22mm male
<b>PEEP Connection</b> .....	30mm male
<b>O<sub>2</sub> Tubing Size</b> .....	6.5ft (198cm) x 6mm OD
<b>O<sub>2</sub> Tubing Inlet Connection</b> .....	15mm x 6mm OD
<b>Reservoir Bag Volume</b> ....	2600ml
<b>Shelf Life</b> .....	5 years

**PerSys Medical**  
5310 Elm St.  
Houston, TX 77081 USA  
1.888.737.7978  
info@ps-med.com

[www.persysmedical.com](http://www.persysmedical.com)

