

# Product Procedure Specifications Sheet

Specification Sheet No	SSN-JAL007-003			
Product Name	Medical Kits			
Product Description	Polyfield Dressing Pack			
Range	Premier			
Manufacturer	Allmed Medical Products Co., Ltd.			
Packaging	80pcs			
Sterilisation	ETO			
Issue Date	28 <sup>th</sup> Feb. 2025	Version No	R02	

Intended use	Intended to be used to assist in minor woundcare procedures
--------------	---

7075AF Pack content				
Item	Product code	Quantity	Dimensions	Material
Non-Woven Swabs	B02213-0275	10	10 X 10CM 4 PLY	Polyester and viscose
Hand Towel	Y22093	1	38 X 40CM	Fiber Paper
White Waste Bag with Tie Handles	P441162	1	300x520mm	PE
Tissue Backed Drape	Y22031	2	50 X 50CM	Polyethylene and fibre paper
Apron	Y25012	1	70 X 135 CM	PE
Measuring Tape	Y33070	1	20CM	Paper
1 Part Tray	Y21113	1	103x117x19mm	PP
Nitrile Gloves No Wallet	Y20027	1 Pair	SMALL size (Palm 80mm)	Nitrile
Tissue Backed Sterile Field	Y22051	1	50 X 50CM	PE+ White paper

7076AF Pack content				
Item	Product code	Quantity	Dimensions	Material
Non-Woven Swabs	B02213-0275	10	10 X 10CM 4 PLY	Polyester and viscose
Hand Towel	Y22093	1	38 X 40CM	Fiber Paper
White Waste Bag with Tie Handles	P441162	1	300x520mm	PE
Tissue Backed Drape	Y22031	2	50 X 50CM	Polyethylene and fibre paper
Apron	Y25012	1	70 X 135 CM	PE
Measuring Tape	Y33070	1	20CM	Paper
1 Part Tray	Y21113	1	103x117x19mm	PP
Nitrile Gloves No Wallet	Y20028	1 Pair	Middle size (Palm 95mm)	Nitrile
Tissue Backed Sterile Field	Y22051	1	50 X 50CM	PE+ White paper

7077AF Pack content				
Item	Product code	Quantity	Dimensions	Material
Non-Woven Swabs	B02213-0275	10	10 X 10CM 4 PLY	Polyester and viscose
Hand Towel	Y22093	1	38 X 40CM	Fiber Paper

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004-03	Procedure Pack Specification Sheet	Jack Wong	Racheal Begg	R02	03 <sup>rd</sup> April 2024

White Waste Bag with Tie Handles	P441162	1	300x520mm	PE
Tissue Backed Drape	Y22031	2	50 X 50CM	Polyethylene and fibre paper
Apron	Y25012	1	70 X 135 CM	PE
Measuring Tape	Y33070	1	20CM	Paper
1 Part Tray	Y21113	1	103x117x19mm	PP
Nitrile Gloves No Wallet	Y20026	1 Pair	LARGE size (Palm 110mm)	Nitrile
Tissue Backed Sterile Field	Y22051	1	50 X 50CM	PE+ White paper

<b>Quality Certification</b>	ISO 13485		
<b>Device Classification</b>	Class Is	<b>Rule</b>	4
<b>Conformity Route</b>	Annex II+III and Annex XI-Part A of Regulation (EU) 2017/745		
<b>Storage conditions</b>	Store in a cool dry place away from heat and direct sunlight		
<b>Shelf life</b>	3 Years - under recommended storage conditions		
<b>Contraindications/Warnings</b>	No Known		
<b>Markings</b>	Product Name, REF, Premier Logo, Lot Number, Manufacturing Date, Expiry Date, Medical Device, Single Use, EO Sterile, Single sterile barrier, do not use if package damaged, Keep dry and away from sunlight, CE, ECREP, Manufacturer name and address, Distributor name and address		
<b>Additional information</b>	Products are Latex and phthalate free Products contain no hazardous substances Products contain no medicinal substances		
<i>Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.</i>			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: <a href="https://mediq.co.uk/">https://mediq.co.uk/</a>	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd.  
Registered in England No. 62537

## Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	12/06/2023	1 <sup>st</sup> version	Jack Wong
R02	28/02/2025	2 <sup>nd</sup> version	Elizabeth Ogunkolade

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004-03	Procedure Pack Specification Sheet	Jack Wong	Racheal Begg	R02	03 <sup>rd</sup> April 2024