



Welch Allyn

Pocket Plus LED

Otoscopes,
Ophthalmoscopes and Sets

Maximized performance and durability in a compact set.



Blueberry



Plum



Onyx



Vanilla

Examine your patients in style, with four different colors to choose from, including colorful, protective accessory kit.

WelchAllyn



Develop and perfect your physical assessment skills while making your mark on the medical field.

Pocket Plus LED Otoscope

- The most durable pocket-style otoscope on the market today, based on rigorous drop test versus leading competitors
- Features up to 6 lumens for increased light intensity
- Fiber optics provide cool light with no reflections or obstructions
- Lightest-weight Welch Allyn otoscope, for the comfort and control you need

Pocket Plus LED Ophthalmoscope

- All the apertures and filters as in our full-size ophthalmoscopes, so you don't have to sacrifice any functionality
- Features a neutral density light filter switch, which helps to virtually eliminate corneal glare
- Lightest-weight Welch Allyn ophthalmoscope for portability and ease of use



The Welch Allyn LED Difference

We've optimized Welch Allyn LED lamps for EENT instruments, providing the benefits of LED light, to offer true tissue color, clarity and definitions of the ear canal and fundus

- Rated for up to 20,000 hours of use
- No lamp replacement required

92880-BLU Pocket Plus LED Set, Blueberry with Hard Case
92880-PUR Pocket Plus LED Set, Plum with Hard Case

92880-BLK Pocket Plus LED Set, Onyx with Hard Case
92880-WHT Pocket Plus LED Set, Vanilla with Hard Case



WelchAllyn®

Explore the entire Welch Allyn physical exam line at welchallyn.com/exam or call us at 800.535.6663.

Welch Allyn, Inc.
4341 State Street Road, Skaneateles Falls, NY 13153 USA
(p) 800.535.6663 (f) 315.685.3361

WWW.WELCHALLYN.COM

© 2015 WELCH ALLYN MC12174