




FREESTANDING PHARMACY FRIDGE SOLID DOOR 353L WHITE



- ✓ **+2°C to +8°C Operating Temperature Range**
Conforms to MHRA guidelines
- ✓ **NEW Advanced Control**
Temperature display accurate to 0.1°C
- ✓ **NEW Unique Anti Microbial Handle**
Helps prevent the growth of bacteria and minimise the potential for cross contamination
- ✓ **Fitted Lock as Standard**
Maintains security and ensures the unit can only be accessed by approved personnel
- ✓ **Interior Light**
- ✓ **Fitted Rollers**
For ease of movement
- ✓ **Wire Shelves and Basket**
Enables optimum air flow inside the unit
- ✓ **Audible Warning Alarms**
Alerts user should internal temperature deviate from set parameters



www.facebook.com/Lec-Medical
www.twitter.com/Lec_Medical

SKU CODE: PRODUCT CODE:

444442729 PSR353UK

DESCRIPTION:

Freestanding Pharmacy Fridge Solid Door 353L White



FULL SPECIFICATION – PSR353UK FREESTANDING PHARMACY FRIDGE SOLID DOOR 353L WHITE

Description	Freestanding
SKU Code	444442729
GDPA Code	PSR353UK
Gross Capacity (litres)	353
Net Capacity (litres)	349
Unit External Width (mm)	595
Unit External Depth (mm)	640 (w/o handle), 670 (with handle)
Unit External Height (mm)	1860
Product Weight (kg)	69
Packed Width (mm)	662
Packed Depth (mm)	732
Packed Height (mm)	1932
Packed Weight (kg)	75
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	480 x 445 x 1682 (step: 175d x 265h)
Temperature Range	+2°C to +8°C
Colour	White/Green
External Finish	Painted Steel
Internal Finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS Mix)
Door Finish	Painted Steel
Door Type	Solid
Door Style	Hinged, Reversible
Door Open Depth (mm)	1185
Replacement Door Seal	No
Open Door Sensor	Yes
Lock	Yes, 2 Keys
Adjustable Feet (number & position)	2 Adjustable Front Feet, 2 Rear Rollers
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Number of Shelves	6
Number of Shelf Positions	14
Shelf Size (mm) (WxDxH)	6: 385 x 500 x 6
Shelf Material	Wire
Number of Baskets	1
Basket Material	Wire
Interior Light	Yes
Fan	Yes
Spark Free	No
Auto or Manual Defrost	Automatic
Controller type	Digital, External
Climate Class	SN
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	Yes
HCFC Free	Yes
Ammonia Free	Yes
UK 3 Pin Plug	Yes
Other Plug Types Available	Yes
Voltage Supply	220-240v/50hz
Aux Output	Yes
Insulation Thickness (mm)	60
Insulation Type	Cyclopentane
Anti Microbial Handle	Yes
Optional Drawer System Available	Yes