

Product Specification



Version: 04:00

Date of issue: 27/02/2017

Date of revision: 06/11/2020

Medipal Chlorine Disinfectant Wipes

Medipal Chlorine Disinfectant Wipes are intended for disinfection of hard surfaces of non-invasive medical devices within healthcare and care environments. For professional use only.

Medipal Chlorine Disinfectant Wipes are suitable for a high level of disinfection, proven to kill most microorganisms such as bacteria, spores, yeast and viruses. The wipes have been independently tested to a range of European testing standards (EN) which includes EN16615.



Product Description

Product Name	50 Medipal Chlorine Disinfectant Wipes
Product Code	S570110MPCE
Languages	EN, FR, ES, AR

Pack Type & Dimensions

Count	50
Dimensions	200 x 200mm
Weight	0.710kg
Colour of wipes	White
Wipes material description	50gsm 100% PET Spunlace

Solution

Solution Description	Medipal Chlorine Solution
Active Ingredients Content	Sodium Hypochlorite 5 g/kg
Hazard identification	For further information please refer to the material safety data sheet

Primary Packaging

Pack Type	Flow wrap pack (film laminate)
Lid Description	PP Rigid Lens
Label	Not applicable
Barcode	5025254 028728

Product Specification



Version: 04:00

Date of issue: 27/02/2017

Date of revision: 06/11/2020

Medipal Chlorine Disinfectant Wipes

Secondary Packaging

Case Description	Plain outer case
Barcode	5025254 022399

Packaging & Logistical Information

Product Name	50 Medipal Chlorine Disinfectant Wipes
Product Code	S570110MPCE
Case Quantity	12 Packs/case
Case Dimensions [OD]	360 x 250 x 190
Case Weight	9.70kg
Case/Layer	11
Layers/Pallet	5
Cases/Pallet	55
Pallet Weight	515kg
Storage Temperature	5-25°C
Hazardous for Transport	Not applicable
Commodity code for export	3808942090