

Product Specification

Version: 01:00

Date of issue: 27/08/2025

Date of revision: 27/08/2025



Medipal Chlorine Disinfectant Wipes

Medipal Chlorine Disinfectant Wipes are intended for disinfection of surfaces within healthcare and care environments. For professional use only.

Medipal Chlorine Disinfectant Wipes are suitable for a high level of disinfection, proven to kill most microorganisms such as bacteria, spores, yeast and viruses. The wipes have been independently tested to a range of European testing standards (EN) which includes EN16615.

Product Description

Product Name	5 Large Medipal Chlorine Disinfectant Wipes
Product Code	S573110MP
Languages	EN, FR, PT, ES, NL

Pack Type & Dimensions

Count	5
Dimensions	150 x 200mm
Weight	0.26kg
Colour of wipes	White
Wipes material description	50gsm 100% PET Spunlace

Solution

Solution Description	Medipal Chlorine Solution
Active Ingredients Content	Sodium Hypochlorite 5 g/kg
Hazard identification	For further information please refer to the material safety data sheet

Primary Packaging

Pack Type	Flow wrap pack (film laminate)
Lid Description	Rigid Lens
Label	Not applicable
Barcode	5025254 001097

Product Specification

Version: 01:00

Date of issue: 27/08/2025

Date of revision: 27/08/2025



Medipal Chlorine Disinfectant Wipes

Secondary Packaging

Case Description	Plain outer case
Barcode	5025254 047521

Packaging & Logistical Information

Product Name	5 Large Medipal Chlorine Disinfectant Wipes
Product Code	S573110MP
Case Quantity	24 Packs/case
Case Dimensions [OD]	240 x 170 x 275
Case Weight	5.60kg
Case/Layer	30
Layers/Pallet	5
Cases/Pallet	150
Pallet Weight	411kg
Storage Temperature	5-25°C
Hazardous for Transport	Not applicable
Commodity code for export	3808942090