

# Sharpak 35 PLUS, 60 PLUS, 120 PLUS and 220 PLUS sharps containers

The Sharpak range of blue-top sharps containers are recommended for the segregation, collection and disposal of partially discharged sharps waste which are NOT contaminated with cytotoxic or cytostatic medicines.

This comprehensive range of 10 containers is available from 0.5L through to 26.3L and is accompanied by a range of brackets, trays and stands.

Sharpak sharps containers are approved to ISO 23907:2012 and tested and certified in accordance with the provisions of the United Nations recommendations on the Transportation of Dangerous Goods, Chapter 9, UN3291 Clinical Waste, Unspecified, NOS.

Sharpak sharps containers MUST only be disposed of by incineration.

Please note:- For disposal of sharps contaminated with cytotoxic or cytostatic medicines please refer to the data sheet for the Code PURPLE range of Sharpak HYDRI sharps containers.

MANUFACTURED, TESTED AND CERTIFIED TO ISO 23907:2012

UN TYPE TESTED

BLUE COLOUR CODED CLOSURES TO MEET THE RECOMMENDATIONS ON WASTE SEGREGATION AS DETAILED IN THE DEPARTMENT OF HEALTH'S TECHNICAL MEMORANDUM HTM-07-01

SECURE CLOSURE WITH LOCKING FACILITY

VISIOLOCK INDICATOR FITTED ON THE SHARPAK 35 PLUS, 60 PLUS, 120 PLUS AND 220 PLUS

DISTINCTIVE COLOURING AND CLEAR LABELLING

SEQUENTIAL NUMBERING FOR TRACEABILITY

TOUGH AND DURABLE

QUICK AND EASY TO ASSEMBLE

COST EFFECTIVE

PURPOSE DESIGNED WALL/TROLLEY BRACKETS, KICK TROLLEY, PROCEDURE TRAY AND WALL UNIT

INSULIN PEN NEEDLE REMOVING FEATURE ON SHARPAK 35 PLUS



Patented and trademarked by Nolato.



3.5 LITRES

6 LITRES

12 LITRES


22 LITRES



PRODUCT	MAX CAPACITY	CAPACITY TO FILL LINE	PRODUCT CODE	QTY	HEIGHT	WIDTH	DEPTH
Sharpak 35 PLUS	3.5 litres	2.63 litres	SU03500	40	195 mm	177 mm	206 mm
Sharpak 60 PLUS	6 litres	4.5 litres	SU02600	25	229 mm	220 mm	234 mm
Sharpak 120 PLUS	12 litres	9.6 litres	SU02200	15	325 mm	258 mm	268 mm
Sharpak 220 PLUS	22 litres	17.6 litres	SU02400	8	383 mm	315 mm	327 mm

ACCESSORIES	CODE	QTY	ACCESSORIES	CODE	QTY
Kick Trolley for 120 or 220 PLUS	SM16046	1	Wall & trolley bracket for 35 - 220 PLUS	SM05000	1
Shaped procedure tray for 35 PLUS	SU03501	5	SAFE & SOUND wall unit for 60 PLUS	SS00060	1



For details regarding Sharpak distributors and worldwide availability, please contact:  
 Nolato Jaycare Ltd Walton Road Portsmouth PO6 1TS  
 Tel +44 2392 314711 Fax +44 2392 380314 www.sharpaks.com