



Be sure with us

Pharmacy refrigeration range

www.shoreline-medical.co.uk

Pharmacy range features

- ✓ *Meets RPSGB, NHS and RCVS criteria*
-  *Suitable for storage of pharmaceuticals and vaccines*
-  *Temperature range +2°C to +8°C*
-  *Forced-air cooling for even temperature distribution*
-  *Audio/visual alarm on temperature deviation*
-  *Max/Min temperature memory*
-  *Lockable doors*
-  *Automatic defrosting*
-  *CFC and Ammonia free*
-  *Adjustable shelves*



Shoreline (UK) Ltd reserves the right to change specifications without notice in line with our policy of continual improvement.



Glass door models



Model SM544G

544 litre capacity
5 shelves
Dimensions: 750 x 750 x 1700mm



Model SM364G

335 litre capacity
5 shelves
Dimensions: 600 x 620 x 1640mm



Model SM264G

236 litre capacity
4 shelves
Dimensions: 600 x 620 x 1125mm



Model SM161G

141 litre capacity
3 shelves
Dimensions: 600 x 615 x 835mm



Model SM60G

60 litre capacity
2 shelves
Dimensions: 470 x 455 x 660mm

Dimensions: width x depth x height

Our promise



Solid door models



Model SM544
544 litre capacity
5 shelves
Dimensions: 750 x 750 x 1700mm



Model SM364
335 litre capacity
5 shelves
Dimensions: 600 x 620 x 1640mm



Model SM264
236 litre capacity
4 shelves
Dimensions: 600 x 620 x 1125mm



Model SM1502
151 litre capacity
3 shelves
Dimensions: 550 x 610 x 850mm



Model SM161
141 litre capacity
3 shelves
Dimensions: 600 x 615 x 835mm



Model SM1402
141 litre capacity
3 shelves
Dimensions: 500 x 610 x 850mm



Model SM75
75 litre capacity
2 shelves
Dimensions: 450 x 430 x 790mm



Model SM45
45 litre capacity
2 shelves
Dimensions: 380 x 375 x 715mm

Product features



Castors on large models



User-friendly temperature control



Door security



Adjustable shelf positions



Forced air cooling

Pharmacy refrigeration

* including bottom rear step

[illegible]

Portable pharmacy refrigerators



Model SMP26 Portable

16kgs
33 watts
26 litre capacity
Dimensions: 345 x 510 x 385mm



Model SMP41 Portable

22kgs
33 watts
41 litre capacity
Dimensions: 610 x 385 x 400mm



Model SMP65 Portable

24kgs
42 watts
65 litre capacity
Dimensions: 725 X 510 x 450mm

Dimensions: width x depth x height

Portable pharmacy range features



12 Volt DC and
100-240 AC operation



Meets RPSGB, NHS
and RCVS criteria



Suitable for storage
of pharmaceuticals
and vaccines



Temperature range
+2°C to +8°C



Forced-air cooling
for even temperature
distribution



Audio/visual alarm on
temperature deviation



Max/Min temperature
memory



Automatic defrosting



Finished in durable
moulded exterior

Technical Specifications

Pharmacy refrigeration

* including bottom rear step

	SMP26	SMP41	SMP65
External Dimensions (W x D x Hmm)	345 x 510 x 385	610 x 385 x 400	725 x 510 x 450
Internal Dimensions (W x D x Hmm)	275 x 435 x 255*	525 x 290 x 260*	610 x 360 x 325*
Product Weight (Kgs)	16	18	24
Gross Capacity (Litres)	26	41	65
Net Capacity (Litres)	24	39	63
Temperature Range	+2°C to +8°C	+2°C to +8°C	+2°C to +8°C
Door Style	Hinged Lid	Hinged Lid	Hinged Lid
Replaceable Door Seal	Yes	Yes	Yes
Lock (available on request)	Upon request	Upon request	Upon request
Carry handles	Yes	Yes	Yes
Temperature Display	Yes	Yes	Yes
Temperature Alarm	Yes	Yes	Yes
MAX / MIN Temperature Memory	Yes	Yes	Yes
Shelves	Basket	Basket	Basket
Interior Light	Yes	Yes	Yes
Automatic Defrost	Yes	Yes	Yes
Interior Fan	Yes	Yes	Yes
Climate Class	SN	SN	SN
Refrigerant Type	R134a	R134a	R134a
UK 3-pin plug	Yes & Car Socket	Yes & Car Socket	Yes & Car Socket
Voltage	12v and 100-240v	12v and 100-240v	12v and 100-240v
Frequency	50/60	50/60	50/60
Calibration Certificate	Yes	Yes	Yes
Warranty (UK only)	2 years	2 years	2 years

Only Shoreline offers a free 5 year calibration service and warranty



Because we have been specialists in medical refrigeration for the past 30 years, we know what's important to healthcare professionals, and assurance of reliability and performance has to be essential!

Whether you are storing vaccines or temperate sensitive pharmaceuticals, like retrovirals and reagents, or even biological samples, you have to know that your patients' interests are always protected and expensive drugs will not be damaged or lost.

So we have designed and built for your needs and include a warranty offer and an after-sales service that no other supplier has matched.

Our product range incorporates our specialist 'know-how' in key aspects like temperature control and monitoring, energy efficiency through thick wall insulation and LED lighting, ease of use and product

security, whilst aiming to deliver best value. We have also added helpful features like a min/max temperature memory, automatic defrosting, adjustable shelves, even castors on large models.

And, of course, our range is complaint with all the relevant and international standards including those of NHS, RPSGB and RCVS.

Shoreline (UK) Ltd is an ISO9001 registered company.

For information, technical support and advice, prices and availability, just call us on 01903 733 877 or email: ask@shoreline-medical.co.uk



+44 (0) 1903 733 877 | ask@shoreline-medical.co.uk | www.shoreline-medical.co.uk

Shoreline (UK) Ltd, 16-17 Martello Enterprise Centre, Courtwick Lane, Littlehampton, West Sussex BN17 7PA