


Product Specifications Sheet

Specification Sheet No	SSN-FT1357-001		
Product Name	Wound Care		
Product Description	Wound Closure Strips (365 Strips)		
Range	365 Healthcare		
Manufacturer	Shanghai ISO Medical Products Co Ltd (SIMP)		
Material	Non-woven fabric white 20gsm, Polyester line high tensile, Solvent adhesive 30gsm, Release paper 120gsm		
Sterilisation	EO Sterilisation		
Issue Date	21st July 2025	Version No	R03



Intended use	A single use set of wound closure strips intended to be used for skin closures including superficial lacerations and skin tears.
---------------------	--

Packaging									
Product code	Colour	A (±1) mm	B (±3) mm	C (±1) mm	D (±2) mm	Strips / pouches	Items Per Pack	Qty- Inner	Qty – Shipper Carton
36519064	White	76	36	3	12	5	16	50	800
36519065		76	39	6		3			
36519066		38	39	6		6			
36519067		102	51	6		10			
36519068		102	51	13		6			
36519069		102	57	25		4			

Quality Certification	ISO 13485:2016		
Device Classification	Class IIa	Rule	6
Conformity Route	Medical Devices Directive (93/42/EEC)		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 Years if stored as recommended		
Contraindications/Warnings	None		
Markings	365 Healthcare, Ref, Product name, Size, Lot, Number of strips, Quantity, Barcodes, EC-REP, EO Sterile, Manufacturing date, Expiry date, CE Mark, Single use, Latex free, Medical device, Do not use if package is damaged, keep dry, keep away from sunlight, Manufacturer Name & Address		
Additional information	Latex free and sterile. Do not use if the packaging is broken.		
<i>Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.</i>			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024



365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	04/10/2023	New SSN for Supplier CE/365 Branding	E.Ogunkolade
R02	21/05/2024	Update Mediq Healthcare UK address	Jack Wong
R03	21/07/2025	Update the classification from Class Is to Clas IIa	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024