

## TANGO

### Integrated pediatric and adult spine boards



Tango spine boards are devices designed to lift and immobilize patients with suspected spinal injuries.

#### Specific features

The integration of adult and pediatric spinal boards, offers unmatched versatility

Pediatric spinal board with differentiated occipital areas for the best adaptation to the anatomical characteristics of pediatric patients

The pediatric spinal board can be easily extracted from the adult one, ensuring fast rescue operations and minimizing storage space

The head immobilizer can easily placed thanks to strap bands

#### Technical data

Spinal board dimensions	1830 x 445 x 55 mm
Number of handles	14
Pediatric spinal board weight	3 kg
Adult spinal board weight	8 kg
Weight of complete spinal board	11 kg
Tango load capacity	150 kg
Baby Go load capacity	30 kg
Pediatric spinal board dimensions	1190 x 320 x 45 mm
Pediatric spinal board number of handles	10
Materials	Polyethylene

#### Accessories

ST02102B	Spine Pack – Carry bag
ST02106A	MARK UP – Permanent personalization system
ST02105A	Fix Board – 10 G Wall support for spinal board

Class I MD compliant with UE Reg. 2017/745

TANGO ADULT AND PAEDIATRIC SPINE BOARD



**ST02140C**

**CND Classification** V0804

**Registration number** 101670

Rev.0 (14/06/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)