

STX 519

Adjustable lifting bridles system



These harnesses are accessories for basket stretchers, to be used solely for lifting from a floor-mounted position or for dragging.

Specific features

4 hooks with automatic locking secure
Easy to adjust system

Technical data

Belt length	Maximum 1680 mm Minimum 1250 mm
Belt width	50 / 20 mm
Weight	1.8 ± 0.1 kg
Material	Nylon, aluminium, steel
Fastening type	Spring catch
Number of fastening bands	4
Compatibility	All basket stretchers manufactured by Spencer
Maximum applicable static safety load	450 kg (1000 lbs)

Class I MD compliant with UE Reg. 2017/745

STX 519 ADJUST. LIFTING BRID.
- LIFTING BRIDLE SYS

ST04519C

CND Classification V0899

Registration number 188142



Rev.0 (27/11/2020)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it