

# SPENCER SX-R

## Pickup stretcher



Spencer SX is a pickup stretcher with three blades for lifting and loading of patients on immobilizing or transport devices such as spinal boards, vacuum mattresses or transport stretchers.

### Specific characteristics

- Extremely light and resistant structure
- The possibility to adjust the length confers high versatility
- Foldable to save space when stored
- Smooth surfaces easy to clean
- The numerous handles make the transport easier
- Belts with transport bag

### Standard accessories

- ST00593A K belt 3—three one piece belts STX597

### Optional accessories

- ST05009A SX Bag – Orange transport bag
- ST05018C SXV – 20G vertical fixing system
- ST05019C SXO – 20G horizontal fixing system

### Compliant with

UNI EN 1865-1

### Technical data <sup>1</sup>

|                       |              |
|-----------------------|--------------|
| Shortened length      | 1660 mm      |
| Intermediate length 1 | 1780 mm      |
| Intermediate length 2 | 1900 mm      |
| Intermediate length 3 | 2050 mm      |
| Folded length         | 1190 mm      |
| Thickness folded      | 70 mm        |
| Thickness opened      | 65 mm        |
| Width                 | 440 mm       |
| Material              | Aluminum     |
| Loading capacity      | 170 kg       |
| Weight                | 7,9 ± 0,4 kg |

<sup>1</sup> Length measures are subjected to tolerances of ± 10 mm

Class I Med.Device according to Dir. 93/42/CEE

SPENCER SXR PICK UP STRETCHER  
3 BLADES, W/3 BELTS

ST05001B

CND Classification V9099

Registration number 91816



Rev.3 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)