

ONE STEP. ONE MINUTE. ONE SOLUTION. FOR ALCOHOL-FREE DISINFECTION



Why Choose Sani-Cloth® AF³?

Sani-Cloth® AF³ wipes are alcohol free, ready-to-use low linting, disposable 2 in 1 cleaning and disinfectant wipes for non-invasive medical devices; non-porous hard surfaces. Sani-Cloth® AF³ can be used to pre-clean or decontaminate medical devices prior to sterilisation.



Test Method	Microorganism	Total Contact Time
EN 16615 EN 13727	Bactericidal	1 min (d)
EN 16615 EN 13624	Yeasticidal	1 min (d)
EN 14476	Enveloped Viruses	1 min (d)
EN 14476	Norovirus, Adenovirus	1 min (d)
EN 14348	TB	1 min (d)



Equipment Protection: The largest alcohol-free compatibility database in the world, tested with more than 450 medical devices.



Efficient: 1 minute total contact time including TB & norovirus in dirty conditions.



Skin Integrity: Dermatologically tested. Formulated without PHMB, a known skin sensitiser, to minimise the risk of allergic reactions. Can be used without gloves.

(d) Proven to clean and disinfect on visibly soiled surfaces – no pre-cleaning needed.



Minimising our carbon footprint with local manufacturing at our UK plant.



Sani-Cloth® AF³ is compatible with our Point of Use dispensers for improved compliance!

Product Description	Product Code	NHS CODE	Wipe Size	Classification	PAO	Shelf life
200ct Soft Pack	EAF3001	VJT592	267mm x 210mm	MD Class IIa	6 months	24 months
225ct Bucket	EAF3005	VJT513	260mm x 270mm	MD Class IIa	6 months	36 months
200ct Canister	EAF3002	VJT592	195mm X 200mm	MD Class IIa	6 months	18 months



To find out more write to us at: contact@pdi-intl.com

- contact@pdi-intl.com
- +44 (0) 8081 697 945
- pdi-intl.com
- PDI International