

# PROVEN. POWERFUL. SPORICIDAL.

## COMPLIANT WITH EN 17846



# Why Choose Sani-Cloth® Chlor?

Sani-Cloth® Chlor is a sporicidal, high-level disinfection, disposable cleaning and disinfection wipe for non-invasive medical devices; non porous hard surfaces. These ready to use, low linting wipes can be used to clean and disinfect in one product. Sani-Cloth® Chlor can be used to pre-clean or decontaminate medical devices prior to sterilisation.



ECH1400



Test Method	Microorganism	Total Contact Time
EN 16615 EN 13727	Bactericidal	2 min (d) 30 sec (c)
EN 16615 EN13624	Yeasticidal	2 min (d) 1 min (c)
EN 14476	Enveloped Viruses	1 min (d)
EN 14476	Non-enveloped Viruses	1 min (d)
EN 17846	C.Difficile	1 min (d)



**Equipment Protection:** The largest sporicidal compatibility database in the world, tested with more than 900 medical devices.



**Efficient:** 2 minute total contact time including *C.Difficile* in 1 minute in clean and dirty conditions.



**Free From:** Formulated without PHMB, a known skin sensitiser, to minimise the risk of allergic reactions.

(c) Effective disinfection on visibly clean surfaces. For soiled areas, pre-cleaning is recommended to ensure full efficacy.  
(d) Proven to clean and disinfect on visibly soiled surfaces - no pre-cleaning needed.



Minimising our carbon footprint with European manufacturing at our UK & Italian plants.

Sani-Cloth® Chlor is compatible with our Point of Use dispensers for improved compliance!



A002GRY

Product Description	Product Code	Wipe Size	Classification	PAO	Shelf life
50ct Canister	ECH1400	200mm x 195mm	MD Class IIa	1 month	24 months