

Unifix 3mL

Luer slip syringe



Sterile 3mL Luer slip syringe for single use.

Supplied in boxes of 100 syringes. An outer carton contains 20 boxes, totalling 2,000 syringes (see page 3 for details).

Product code

UX3

Materials

Syringe body: Medical grade Polypropylene (PP 5250T)

Plunger gasket: Isoprene rubber (latex-free)

Syringe barrel and gasket lubricating oil: Medical grade Silicone oil (DC360)

None of these materials contain recycled content

Sterilisation

Sterilised by: Ethylene oxide gas

Conforming to: ISO:11135 Sterilisation of medical devices
(validation and routine control of ethylene oxide sterilisation)

Time to expiry: 5 years from manufacture

Product conformance

Class IIa Medical Device

Fully compliant with ISO:7886-1: Sterile hypodermic syringes for single use

European Medical Device Directive 93/42/EEC (CE Class:1)

Manufacturing compliant with Quality Assurance Standards ISO:13485



Graduations in
0.1mL increments

Compatibility

All Luer fit needles

Sterifilt compatible

Manufacture and distribution

Unifix syringes are manufactured exclusively for Exchange Supplies, UK.

1, Great Western Industrial Centre, Dorchester Dorset DT1 1RD. www.exchangesupplies.org

Manufactured for	Shanghai Kindly Enterprise Development Group Co., Ltd.
Exchange Supplies by:	No. 658 Gaochao Road, 201803 Shanghai, People's Republic of China.

EC representative:	Shanghai International Holding Corp. GmbH (Europe). Eiffestraße 80, 20537 Hamburg, Germany. T: 0049-40-2513175. F: 0049-40-255726
---------------------------	--

KDL UKCA representative:	KDL – Overseas Business Development Director. Unit 14 Cowley Mill Road, Uxbridge, Greater London UB8 2DB.
---------------------------------	--

Packaging

Item	Weight (approx.)	Material(s)
Blister	4 g	Printed surface: Medical grade 63gsm coated paper Blister: 0.09mm thickness medical complex film
Box 100	463g	E corrugated card
Carton 2,000	12.95 kg	Uncoated 630gsm AB flute double wall card. Consisting of: 150gsm 3A Kraft Test Liner printed outer layer, with 120gsm and 100gsm A grade Senior corrugated and 140gsm B grade Kraft Test inner Liner. Recyclable.

None of these materials contain recycled content.



Blister
138.6mm x 39.5mm
GTIN: 06930197851137



Box 100
210mm x 135mm x 140mm
GTIN: 16930197851134



Carton 2,000
710mm x 300mm x 460mm
GTIN: 26930197851131

EXCHANGE
S U P P L I E S

1 Great Western Industrial Centre, Dorchester,
Dorset DT1 1RD | T: +44 (0) 1305 262244