


Technical **information**



Designed with safety in mind

The new **Unifine® Syringes** feature clear graduations and smooth plunger action to promote accurate dosage, easy operation, and minimal user error to support healthcare professionals and patient care in a range of medical scenarios.



**Available in three tip options
and multiple barrel sizes**



Device specifications

Size	1mL	2mL	3mL	5mL	10mL	20mL	30mL	50-60mL
Graduation	0.1mL	0.5mL	0.5mL	1mL	1mL	5mL	10mL	10mL
Luer lock central tip								
Device dimensions	98.90 x 20.58mm	87.58 x 22.86mm	87.58 x 22.86mm	85.66 x 26.62mm	100.80 x 30.10mm	127.88 x 37.90mm	136.70 x 41.36mm	145.30 x 47.16mm
Device weight	2.74g	3.46g	3.46g	4.61g	7.16g	12.45g	17.77g	28.79g
Luer slip central tip								
Device dimensions	97.84 x 19.12mm	89.16 x 22.84mm	89.16 x 22.84mm	86.52 x 26.60mm	102.12 x 29.90mm	–	–	–
Device weight	2.48g	3.32g	3.32g	4.37g	7.09g	–	–	–
Luer slip side tip								
Device dimensions	–	–	–	–	113.90 x 30.80mm	129.12 x 38.38mm	136.94 x 38.70mm	160.00 x 52.56mm
Device weight	–	–	–	–	7.96g	13.27g	16.66g	34.75g
Catheter tip								
Device dimensions	–	–	–	–	–	–	–	170.00 x 46.90mm
Device weight	–	–	–	–	–	–	–	29.72g
Device material	Syringe barrel				Polypropylene			
	Plunger							
	Plunger tip				Latex-free rubber			
Sterilisation method	Ethylene Oxide (EtO)							
Regulatory	CE mark, ISO13485, FDA 510(k)							
Shelf life	5 years							
Intended use	The Unifine® Syringe is intended for general purpose injection and aspiration							

Packing and shipping configuration

Boxes

The syringes are packaged into boxes and then packed into cases, which are then palletised for transportation.

Size	1mL	2mL	3mL	5mL	10mL	20mL	30mL	50-60mL
Devices per box	100	100	100	100	100	50	40	50/25

Luer lock central tip

Box size (L x W x H)	20.5 x 14 x 9.5cm	21 x 14 x 14cm	21 x 14 x 14cm	21 x 13 x 20cm	27 x 15.5 x 21cm	22 x 18 x 21.5cm	23.4 x 14.9 x 18cm	39.2 x 21 x 21cm
Box weight	0.390kg	0.550kg	0.550kg	0.778kg	1.175kg	0.913kg	0.563kg	1.833kg

Luer slip central tip

Box size (L x W x H)	21.5 x 14 x 10cm	21 x 14 x 14cm	21 x 14 x 14cm	21 x 13 x 20cm	27 x 15.5 x 21cm	–	–	–
Box weight	0.383kg	0.539kg	0.539kg	0.756kg	1.175kg	–	–	–

Luer slip side tip

Box size (L x W x H)	–	–	–	–	27 x 15.5 x 21cm	22 x 18 x 21.5cm	23.4 x 14.9 x 18cm	39.2 x 21 x 21cm
Box weight	–	–	–	–	1.175kg	0.788kg	0.563kg	1.833kg

Catheter tip

Box size (L x W x H)	–	–	–	–	–	–	–	23.4 x 14.9 x 18cm
Box weight	–	–	–	–	–	–	–	0.938kg

Cases

Boxes are packed into cardboard cases, which are then palletised for transportation.

Size	1mL	2mL	3mL	5mL	10mL	20mL	30mL	50-60mL
Devices per case	3000	1800	1800	1800	1200	400	640	300/400

Luer lock central tip

Case size (L x W x H)	65.5 x 30 x 52cm	44.5 x 44.5 x 44.5cm	44.5 x 44.5 x 44.5cm	65.5 x 41.5 x 42.5cm	56.5 x 49 x 44.5cm	46.5 x 38.5 x 45.5cm	61.5 x 38 x 48.5cm	65.5 x 41 x 44.5cm
Case weight	12.7kg	10.9kg	10.9kg	15kg	15.1kg	8.3kg	10kg	12kg

Luer slip central tip

Case size (L x W x H)	63 x 29 x 49.5cm	44.5 x 44.5 x 44.5cm	44.5 x 44.5 x 44.5cm	65.5 x 41.5 x 42.5cm	56.5 x 49 x 44.5cm	–	–	–
Case weight	12.5kg	10.7kg	10.7kg	14.6kg	15.1kg	–	–	–

Luer slip side tip

Case size (L x W x H)	–	–	–	–	56.5 x 49 x 44.5cm	46.5 x 38.5 x 45.5cm	61.5 x 38 x 48.5cm	65.5 x 41 x 44.5cm
Case weight	–	–	–	–	15.1kg	7.3kg	10kg	12kg

Catheter tip

Case size (L x W x H)	–	–	–	–	–	–	–	61.5 x 38 x 48.5cm
Case weight	–	–	–	–	–	–	–	16kg

Order codes

		Product	Box qty.	GTIN box	Case qty.	GTIN case
Unifine® Syringes Luer lock central tip	1mL	P005279A1	100	16941845295286	3000 (30 boxes)	66941845295281
	2mL	P005467A1	100	16941845295279	1800 (18 boxes)	66941845295274
	3mL	P005280A1	100	16941845295262	1800 (18 boxes)	66941845295267
	5mL	P005221A1	100	16941845295255	1800 (18 boxes)	66941845295250
	10mL	P005222A1	100	16941845295248	1200 (12 boxes)	66941845295243
	20mL	P005223A1	50	16941845295231	400 (8 boxes)	66941845295236
	30mL	P005224A1	40	16941845295224	640 (16 boxes)	66941845295229
	50-60mL	P005225A1	50	16941845295217	300 (6 boxes)	66941845295212
Unifine® Syringes Luer slip central tip	1mL	P005275A1	100	16941845295385	3000 (30 boxes)	66941845295380
	2mL	P005276A1	100	16941845295378	1800 (18 boxes)	66941845295373
	3mL	P005464A1	100	16941845295361	1800 (18 boxes)	66941845295366
	5mL	P005218A1	100	16941845295354	1800 (18 boxes)	66941845295359
	10mL	P005219A1	100	16941845295347	1200 (12 boxes)	66941845295342
Unifine® Syringes Luer slip side tip	10mL	P005465A1	100	16941845295330	1200 (12 boxes)	66941845295335
	20mL	P005220A1	50	16941845295323	400 (8 boxes)	66941845295328
	30mL	P005277A1	40	16941845295316	640 (16 boxes)	66941845295311
	50-60mL	P005278A1	50	16941845295309	300 (6 boxes)	66941845295304
Unifine® Syringes Catheter tip	50-60mL	P005466A1	25	16941845295293	400 (16 boxes)	66941845295298

Ordering information

Contact details:

Joanne Thompson +44 (0) 1608 645555

United Kingdom, Republic of Ireland, Denmark, Finland, Norway and Sweden.

Thomas Rossbach +49 (0) 6026 97750

Germany, Austria, Switzerland, Poland, Czech Republic, Netherlands, Luxembourg, Belgium, Bulgaria, Croatia, Estonia, Hungary, Latvia, Lithuania, Romania, Slovakia, Slovenia, Iceland and Liechtenstein.

Eric Del Cotto +(33) 2 32 51 88 70

France, Italy, Greece, Malta, Portugal and Spain.

Ordering lead times:

Please contact your local Owen Mumford representative for information.

Pricing:

Please contact your local Owen Mumford representative for information.