

Viatom Checkme O2

*Smart wrist oximeter with memory function
compatible with Visi-Download and with a two year warranty*



The Viatom Checkme O2 is the next generation pulse oximeter for overnight monitoring of SpO2 and pulse rate. Compatible with Stowood's Visi-Download, the pulse oximeter uses a rechargeable battery and uses a revolutionary patented 'ring' pulse oximeter, for the best patient experience we have yet encountered.

- Four night recording memory
- Easy-to-clean, reusable ring sensor or other sensors available
- Rechargeable battery
- Compatible with industry-standard Stowood Visi-Download (purchase may be required)
- Soft, comfortable silicone wrist strap
- Nail varnish? No problem! Sensor fits around base of thumb away from nails
- The Viatom Checkme O2 is the next generation wrist-worn pulse oximeter.
- Designed from the bottom up for the best patient experience, Stowood have partnered with the manufacturer to bring Visi-Download compatibility to this innovative new monitor.
- With a simple wrist strap, single button, and download cable, this device is easy to use and to clean.
- The reusable ring sensor is the device's secret weapon- incredibly comfortable, you will probably forget you are wearing a monitoring device at all.




Small, light, with 24 hours battery and room for four recordings. Compatible with Visi-Download via the download cable (supplied separately).

Patented Ring sensor (included) never comes off the patient's thumb, is comfortable and soft, and should lead to a better night's sleep than previous devices.

Technical specification

Recorded channels	SpO ₂ , pulse rate, motion
Dimensions	44 x 25 x 15 mm
Weight including batteries	12 g (main unit)
Power	Internal rechargeable batteries. 24 hours expected for typical recording parameters
Sensors	Ring sensor included as standard. Soft tip sensor, paediatric soft sensor, single patient foam and neonatal non-woven SpO ₂ sensors available
Data storage	Internal memory holds four nights' recordings
Setup & signal check	Via OLED screen
Software	Visi-Download as used for oximeters & capnometers. Reports available in six languages
CE	Viatom holds CE 0197
Warranty	24 months from Stowood on main unit, 6 months on probes



Shenzhen Viatom Technology Co., Ltd.
 4E, Building 3, Tingwei industrial Park, Baoan District, 518100 Shenzhen, P.R.China

© 2018 Stowood Scientific Instruments Ltd. All rights reserved. All product names appearing on this document are trademarks or registered trademarks owned, licenced to, promoted or distributed by Stowood Scientific Instruments Ltd, its subsidiaries or affiliates

Stowood Scientific Instruments Ltd,
 Common Road, Beckley,
 Oxford OX3 9UP, GB
 Tel +44 (0)1865 358860
 www.stowood.com

