



WIRELESS ALERTA BUTTON

The **Alerta Wireless Nurse Call Button** is supplied with a wall mountable holder that can be mounted to a wall or bedside unit for holding the **Wireless Nurse Call Button** meaning it is always in one place if the user is bed bound, also making it easy for those affected by dementia to remember its location.

The **Alerta Wireless Nurse Call Button** is also supplied with a safety lanyard for use over the neck, and a wrist strap enabling it to be wrist worn. The button can be mounted in the wall holder with the wrist strap still attached.

KEY FEATURES & BENEFITS

- > Can be strapped round the wrist
- > Provided with a lanyard to put over the patients neck
- > Wall mountable
- > Wristband material: Silicon
- > Wireless connect-and-play system
- > No cables, extra versatile
- > 95m wireless range from transmitter to receiver
- > Takes 1x CR2450 battery
- > System includes: 1x Alerta Button (must be used with Alerta Wall Point Receiver - sold separately)

