

# SHARPSGUARD®

## Your Responsibility




Assemble correctly. "Click" the lid on all around firmly and fill in the label details.

Place on bracket or take to point of use in a POUDS® tray.

Ensure container door is in open position before use (as shown below).

Lock the door when contents reach the fill line.

Complete the label and/or tag the container.



**Remember always dispose of sharps at the points of use!**

Dispose of containers in your nearest sharps wheeled bin or designated storage area. Under no circumstances mix sharps containers with clinical waste bags.

If you have questions or need further information, please contact your local Area Manager or our Customer Service Team on +44 1706 754980 or [info.healthcareuk@mauserpackaging.com](mailto:info.healthcareuk@mauserpackaging.com) and visit our website [www.daniels.co.uk](http://www.daniels.co.uk).



REDEFINING SUSTAINABILITY

[www.daniels.co.uk](http://www.daniels.co.uk)



[info.healthcareuk@mauserpackaging.com](mailto:info.healthcareuk@mauserpackaging.com)

Unless otherwise stated, the information contained herein is at the specific request of the user and has been made available by Daniels Healthcare.

Your Responsibility 01.2023