



## Choose the stethoscope that helps you achieve your best.

3M™ Littmann® Stethoscopes offer superb acoustic performance, consistently high quality, and outstanding value. Nearly all Littmann stethoscopes use tuneable technology, a 3M invention that gives you the ability to hear different frequency sounds without turning over and repositioning the chestpiece.

## Experience the power of tuneable technology.



**LIGHT CONTACT:**  
Hold the chestpiece with light pressure to hear low frequency sounds.



**FIRM PRESSURE:**  
Press the chestpiece a little more firmly to hear higher frequency sounds.


Lightweight		Classic				Cardiology			Electronic	
<p><b>6</b></p> <p>Acoustic Performance</p> <p>Ideal for taking blood pressure readings and making limited patient assessments.</p>		<p><b>7-8</b></p> <p>Acoustic Performance</p> <p>Hard-working, dependable, clinical tools for physical assessment and monitoring of your patients.</p>				<p><b>9-10</b></p> <p>Acoustic Performance</p> <p>Designed to meet the needs of demanding clinicians. Provides outstanding acoustic sensitivity to help you perform detailed diagnostic auscultation.</p>			<p><b>10+</b></p> <p>Acoustic Performance</p> <p>High-performance auscultation. Hear more of what you need to hear with amplification, ambient noise reduction, and optional connectivity.</p>	
Lightweight II S.E. Stethoscope	Select Stethoscope	Classic II Paediatric and Infant Stethoscopes	Classic II S.E. Stethoscope	Classic III™ Stethoscope	Master Classic II™ Stethoscope	Cardiology III™ Stethoscope	Cardiology IV™ Stethoscope	Master Cardiology™ Stethoscope	Electronic Stethoscope Model 3100	Electronic Stethoscope Model 3200
<p>Double sided chestpiece 118 grams Tuneable diaphragm 2 year warranty Chestpiece: 2.12" / 5.4cm</p>	<p>Single sided chestpiece 120 grams Tuneable diaphragm 2 year warranty Chestpiece: 1.75" / 4.4cm</p>	<p>Double sided chestpiece 105, 95 grams 3 year warranty Chestpiece: Paediatric 1.31" / 3.3cm Infant 1.06" / 2.7cm 0.75" / 1.9cm</p>	<p>Double sided chestpiece 135 grams Tuneable diaphragm 3 year warranty Chestpiece: 1.75" / 4.4cm</p>	<p>Double sided chestpiece 150 grams Single-piece tuneable diaphragm on each side of the chestpiece Convertible open bell 5 year warranty Chestpiece: Adult side 1.7" / 4.3cm Paediatric side 1.3" / 3.3cm</p>	<p>Single sided chestpiece 160 grams Tuneable diaphragm 3 year warranty Chestpiece: 1.75" / 4.4cm</p>	<p>Double sided chestpiece 180 grams Available in 22" / 56cm and 27" / 69cm Two tuneable diaphragms Dual-lumen tubing 5 year warranty Chestpiece: Adult side 1.7" / 4.3cm Paediatric side 1.3" / 3.3cm</p>	<p>Double sided chestpiece 177 grams Available in 27" / 69cm Single-piece tuneable diaphragm on each side of the chestpiece Dual-lumen tubing 7 year warranty Chestpiece: Adult side 1.7" / 4.3cm Paediatric side 1.3" / 3.3cm</p>	<p>Single sided chestpiece 185 grams Available in 22" / 56cm and 27" / 69cm Tuneable diaphragm Dual-lumen tubing 7 year warranty Chestpiece: 2.0" / 5.1cm</p>	<p>185 grams 2 year warranty Chestpiece: 2" / 5.1cm</p>	<p>Bluetooth® wireless technology Recording and playback 185 grams 2 year warranty Chestpiece: 2" / 5.1cm</p>

www.littmann.co.uk

Bluetooth® is a registered trademark of Bluetooth SIG. 3M, Littmann, the L Littmann logo, L, Master Cardiology, Cardiology III, Cardiology IV, Master Classic II, Classic III, and the shape of the stethoscope chestpiece are trademarks of 3M. Used under license in Canada. Please recycle. Printed in UK. © 3M 2016. All rights reserved. 70-2011-5462-5



# 3M™ Littmann® Stethoscope Ordering Information and Tube Colour Guide



		Black	Burgundy	Navy Blue	Hunter Green	Plum	Chocolate	Grey	Pine Green	Ceil Blue	Purple	Raspberry	Caribbean Blue	Orange	Seafoam Green	Red	Lilac	Royal Blue	Ocean Blue	Peach	Pearl Pink	Coral	Dark Olive Green	Sand	Lavender	Turquoise	Lemon-Lime	Lime Green	Emerald	Rose Pink	
<b>Model 3200</b>	69 cm	3200BK27	3200BU	3200NB				3200T																							
<b>Model 3100</b>	69 cm	3100BK27	3100BU	3100NB																											
<b>Master Cardiology™</b>	69 cm	2160 2161* 2175* 2176*	2163	2164	2165 2183*	2167							2178										2182*								
<b>Cardiology IV™</b>	69 cm	6152	6153	6154	6155	6166						6158	6157																	6161	
<b>Cardiology III™</b>	69 cm	3128 3131BE¹ 3128BRS² 3157SM⁴ 3152RBW⁶	3129	3130	3134	3135	3137 3137CPR²	3136		3146		3148	3138			3140						3149	3166*		3158*	3159*	3161*			3164	
<b>Master Classic II™</b>	69 cm	2144L 2141*	2146	2147	2632				2634	2633	2143		2630											2636*							
<b>Classic III™</b>	69 cm	5620 5803* 5811*	5627	5622	5624	5831	5809³	5621		5630		5626 5806*	5623 5807*	5629							5633		5812*		5832	5835	5839	5829	5840	5639	
<b>Classic II S.E.</b>	71 cm	2201	2211	2205						2813		2210								2822	2817										
<b>Classic II Paediatric</b>	71 cm	2113		2123					2131			2122	2119 2153*	2155		2113R	2136			2154	2115										
<b>Classic II Infant</b>	71 cm	2114										2125 2157*	2124	2179		2114R	2156			2158	2120										
<b>Select</b>	71 cm	2290	2293						2305	2301	2294	2296	2291					2298	2306	2310	2292										
<b>Lightweight II S.E.</b>	71 cm	2450	2451							2454			2452		2455		2453				2456										

Teaching stethoscopes, catalogue numbers 2138 and 2139 as well as veterinary stethoscope, catalogue number 1392 are also available.

Colours shown are approximate.



**3M Health Care  
Littmann® Brand Products**  
3M United Kingdom PLC  
3M Centre, Cain Road, Bracknell,  
RG12 8HT. Registered in England  
and Wales No. 1123045  
www.3m.co.uk/healthcare  
www.littmann.co.uk

**3M United Kingdom PLC**  
Morley Street  
Loughborough  
Leicestershire  
LE11 1EP

3M, Littmann, the L Littmann logo, L, Master Cardiology, Cardiology III, Cardiology IV, Master Classic II, Classic III, and the shape of the stethoscope chestpiece are trademarks of 3M. Used under license in Canada. Please recycle. Printed in UK. © 3M 2016. All rights reserved. 70-2011-5462-5

### Special chestpiece/headset finishes are only available in a limited number of tubing colours.



**1** Black Edition – black finish chestpiece and eartubes.



**2** Brass Edition – brass finish chestpiece and eartubes.



**3** Copper Edition – copper finish chestpiece and eartubes.



**4** Smoke Edition – smoke finish chestpiece and eartubes.



**5** Rainbow Edition – rainbow finish chestpiece and brass finish eartubes.



**6** Rainbow Edition – rainbow finish chestpiece and black finish eartubes.